

Product Information

Product SKU:	CAB14878	Gene ID:	28512	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	22kDa	Conjugate:	Unconjugated
Calculated MW:	22kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable GTPase activating protein binding activity. Predicted to be involved in I-kappaB kinase/NF-kappaB signaling. Predicted to act upstream of or within several processes, including Ral protein signal transduction; lung alveolus development; and surfactant homeostasis. Located in cytosol and endoplasmic reticulum.

**Recommended Dilution:** WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

**Synonyms:** KBRAS1; kappaB-Ras1; NKIRAS1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 118-192 of human NKIRAS1 (NP\_065078.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.