

Product Information

Product SKU:	CAB14879	Gene ID:	28955	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	15kDa	Conjugate:	Unconjugated
Calculated MW:	10kDa	Isotype:	IgG

Immunogen Information

**Background:** DEXI (Dexi Homolog) is a Protein Coding gene. Diseases associated with DEXI include Fibrosclerosis Of Breast and Immunoglobulin A Deficiency 1.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	MYLE; DEXI
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human DEXI (NP_054734.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.