

Product Information

Product SKU:	CAB14885	Gene ID:	30061	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	63kDa	Conjugate:	Unconjugated
Calculated MW:	63kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a cell membrane protein that may be involved in iron export from duodenal epithelial cells. Defects in this gene are a cause of hemochromatosis type 4 (HFE4).

Recommended Dilution:	WB,1:100 - 1:500
Synonyms:	FPN; FPN1; HFE4; MTP1; IREG1; MST079; MSTP079; SLC11A3; SLC40A1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 200-270 of human SLC40A1 (NP_055400.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.