

SLC40A1 Rabbit Polyclonal Antibody



CAB14885

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

62kDa

Applications:

IHC

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

The protein encoded by this gene is a cell membrane protein that may be involved in iron export from duodenal epithelial cells. Defects in this gene are a cause of hemochromatosis type 4 (HFE4).

Immunogen information

Gene ID:

30061

Uniprot

Q9NP59

Synonyms:

SLC40A1; FPN1; HFE4; IREG1; MST079; MSTP079; MTP1; SLC11A3

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 200-300 of human SLC40A1 (NP_055400.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
