

## CAB14887

---

### Product Information

<b>Product SKU:</b>	CAB14887	<b>Gene ID:</b>	51100	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	41kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	41kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a SRC homology 3 domain-containing protein. The encoded protein interacts with the proapoptotic member of the Bcl-2 family, Bcl-2-associated X protein (Bax) and may be involved in regulating apoptotic signaling pathways. This protein may also be involved in maintaining mitochondrial morphology. Alternate splicing results in multiple transcript variants.

**Recommended Dilution:** WB,1:500 - 1:2000 IHC-P,1:50 - 1:200

**Synonyms:** Bif-1; CGI-61; PPP1R70; dJ612B15.2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 120-365 of human Bif-1 (NP\_057093.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.