

DDX49 Rabbit Polyclonal Antibody



CAB14900

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

28kDa/54kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

54555

Uniprot

Q9Y6V7

Synonyms:

DDX49; R27090_2; Dbp8

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 344-483 of human DDX49 (NP_061943.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
