

Product Information

Product SKU:	CAB14926	Gene ID:	79180	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	27kDa	Conjugate:	Unconjugated
Calculated MW:	27kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables cadherin binding activity. Predicted to be located in plasma membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	SWS1; EFHD2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-78 of human EFHD2 (NP_077305.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.