

# EFHD2 Rabbit Polyclonal Antibody



CAB14926

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

27kDa

### Calculated MW:

26kDa

### Applications:

WB IHC

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50  
- 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

79180

#### Uniprot

Q96C19

#### Synonyms:

EFHD2; SWS1

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-78 of human EFHD2 (NP\_077305.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

