

CAB14930

Product Information

Product SKU:	CAB14930	Gene ID:	79625	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	65kDa	Conjugate:	Unconjugated
Calculated MW:	65kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable heparin binding activity. Involved in several processes, including cellular response to hypoxia; negative regulation of apoptotic process; and nitric oxide mediated signal transduction. Located in extracellular space. Implicated in hypogonadotropic hypogonadism.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	HH25; NORD; C4orf31; NDNF
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 20-240 of human NDNF (NP_078850.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.