

Product Information

Product SKU:	CAB14937	Gene ID:	81569	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	41-50kDa	Conjugate:	Unconjugated
Calculated MW:	41kDa	Isotype:	IgG

Immunogen Information

Background: Involved in epithelial cell differentiation. Predicted to be located in cytoplasm. Predicted to be part of dynactin complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CT57; ACTL8
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-366 of human ACTL8 (NP_110439.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.