

DDX59 Rabbit Polyclonal Antibody



CAB14940

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

64kDa/68kDa

Applications:

WB IHC

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

83479

Uniprot

Q5T1V6

Synonyms:

DDX59; OFD5; ZNHIT5


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human DDX59 (NP_001026895.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded rat testis using DDX59 antibody (CAB14940) at dilution of 1:100 (40x lens).