

NTNG2 Rabbit Polyclonal Antibody



CAB14945

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

47kDa/59kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

84628

Uniprot

Q96CW9

Synonyms:

NTNG2; LHLL9381; Lmnt2; NTNG1; bA479K20.1; netrin-G2

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 18-300 of human NTNG2 (NP_115925.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
