
Product Information

Product SKU:	CAB14952	Gene ID:	112858	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	28kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: Enables p53 binding activity and protein serine/threonine kinase activity. Involved in protein phosphorylation. Located in cytoplasm and nucleus. Part of EKC/KEOPS complex. Implicated in Galloway-Mowat syndrome 4.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: PRPK; BUD32; TPRKB; GAMOS4; Nori-2; Nori-2p; dJ101A2; C20orf64; TP53RK

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-253 of human TP53RK (NP_291028.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.