

CAB14956

Product Information

Product SKU:	CAB14956	Gene ID:	115548	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	89kDa	Conjugate:	Unconjugated
Calculated MW:	89kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables identical protein binding activity. Involved in clathrin coat assembly and clathrin-dependent endocytosis. Located in clathrin-coated pit and clathrin-coated vesicle.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	FCHO2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human FCHO2 (NP_620137.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.