

AMICA1 Rabbit Polyclonal Antibody



CAB14958

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

29kDa/40kDa/43kDa/44kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

120425

Uniprot

Q86YT9

Synonyms:

JAML; AMICA; AMICA1; CREA7-1; CREA7-4; Gm638

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 285-384 of human AMICA1 (NP_694938.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
