

Product Information

Product SKU:	CAB14960	Gene ID:	126382	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	15kDa	Conjugate:	Unconjugated
Calculated MW:	16kDa	Isotype:	IgG

Immunogen Information

Background: Located in nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	TRA16; NR2C2AP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-139 of human NR2C2AP (NP_795361.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.