

## CAB14964

---

### Product Information

<b>Product SKU:</b>	CAB14964	<b>Gene ID:</b>	140838	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	23kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	28kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables N-acetylneuraminase-9-phosphatase activity. Involved in N-acetylneuraminase biosynthetic process. Predicted to be located in cytosol.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** HDHD4; C20orf147; dJ694B14.3; NANP

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-248 of human NANP (NP\_689880.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.