

Product Information

| | | | | | |
|--------------|----------|---------------|--------|-------------|-------------|
| Product SKU: | CAB14971 | Gene ID: | 161835 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Mouse |

Additional Information

| | | | |
|----------------|-------|------------|--------------|
| Observed MW: | 70kDa | Conjugate: | Unconjugated |
| Calculated MW: | 66kDa | Isotype: | IgG |

Immunogen Information

Background: FSIP1 (Fibrous Sheath Interacting Protein 1) is a Protein Coding gene. Diseases associated with FSIP1 include Syndromic X-Linked Intellectual Disability Type 10 and Chromosome 3Q29 Duplication Syndrome.

| | |
|-----------------------|---|
| Recommended Dilution: | WB,1:500 - 1:2000 |
| Synonyms: | HSD10; FSIP1 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human FSIP1 (NP_689810.3). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3. |