

# DDX53 Rabbit Polyclonal Antibody



CAB14973

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

### Calculated MW:

71kDa

### Applications:

IHC

### Reactivity:

Mouse

## Antibody Information

### Recommended dilutions:

IHC 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

This intronless gene encodes a protein which contains several domains found in members of the DEAD-box helicase protein family. Other members of this protein family participate in ATP-dependent RNA unwinding.

## Immunogen information

### Gene ID:

168400

### Uniprot

Q86TM3

### Synonyms:

DDX53; CAGE; CT26

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human DDX53 (NP\_874358.2).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---