

Product Information

Product SKU:	CAB14976	Gene ID:	200010	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	74kDa	Conjugate:	Unconjugated
Calculated MW:	74kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable glucose:sodium symporter activity. Predicted to be involved in sodium ion transport. Located in extracellular exosome.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	SGLT4; SLC5A9
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 230-320 of human SLC5A9 (NP_001011547.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.