

ZNF680 Rabbit Polyclonal Antibody



CAB14980

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

14kDa/61kDa

Applications:

IHC

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

340252

Uniprot

Q8NEM1

Synonyms:

ZNF680

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 415-485 of human ZNF680 (NP_848653.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
