

NPBWR1 Rabbit Polyclonal Antibody



CAB15061

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

36kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

2831

Uniprot

P48145

Synonyms:

NPBWR1; GPR7

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 150-250 of human NPBWR1 (NP_005276.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
