
Product Information

Product SKU:	CAB15096	Gene ID:	6231	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	15kDa	Conjugate:	Unconjugated
Calculated MW:	13kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S26E family of ribosomal proteins. Mutations in this gene are found in Diamond-Blackfan anemia 10. There are multiple processed pseudogenes of this gene dispersed through the genome.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: S26; eS26; DBA10; RPS26

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RPS26 (NP_001020.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.