

# KSR1 Rabbit Polyclonal Antibody



CAB15118

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

110kDa

### Calculated MW:

84kDa/87kDa/97kDa/102kDa

### Applications:

WB IF

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:100

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

8844

#### Uniprot

Q8IVT5

#### Synonyms:

KSR1; KSR; RSU2

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 320-420 of human KSR1 (NP\_055053.1).

#### Storage:


Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot analysis of extracts of various cell lines, using KSR1 antibody (CAB15118) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 10s.



Immunofluorescence analysis of C6 cells using KSR1 Polyclonal Antibody (CAB15118) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.