

RABGAP1L Rabbit Polyclonal Antibody



CAB15124

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

8kDa, 29kDa, 41-43kDa,
68kDa, 90-94kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

9910

Uniprot

B7ZAP0/Q5R372

Synonyms:

RABGAP1L; HHL; TBC1D18

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 50-210 of human RABGAP1L (NP_055672.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
