

Product Information

Product SKU:	CAB15154	Gene ID:	29909	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	37kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable G protein-coupled purinergic nucleotide receptor activity. Predicted to be involved in G protein-coupled receptor signaling pathway. Predicted to act upstream of or within negative regulation of myeloid cell differentiation. Predicted to be located in plasma membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	H963; GPR171
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human GPR171 (NP_037440.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.