

CAB15169

Product Information

Product SKU:	CAB15169	Gene ID:	55342	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	74kDa	Conjugate:	Unconjugated
Calculated MW:	74kDa	Isotype:	IgG

Immunogen Information

Background: Enables RNA binding activity. Predicted to be involved in spermatogenesis. Predicted to act upstream of or within mechanosensory behavior and spermatid development. Located in nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	p74; SPNR; ILF3L; HEL162; STRBP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 40-190 of human STRBP (NP_060857.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.