

SHTN1 Rabbit Polyclonal Antibody



CAB15179

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

23kDa, 52- 73kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

57698

Uniprot

A0MZ66

Synonyms:

SHTN1; KIAA1598; shootin-1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SHTN1 (NP_001120683.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
