

AEN Rabbit Polyclonal Antibody



CAB15182

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

36kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

64782

Uniprot

Q8WTP8

Synonyms:

AEN; ISG20L1; pp12744

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human AEN (NP_073604.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
