

Product Information

Product SKU:	CAB15193	Gene ID:	80219	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	27kDa	Conjugate:	Unconjugated
Calculated MW:	27kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable ubiquinone binding activity. Predicted to be involved in cellular respiration and ubiquinone biosynthetic process. Predicted to be located in mitochondrial inner membrane. Predicted to be active in mitochondrion.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	COQ10B
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 38-150 of human COQ10B (NP_079423.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.