

TCTE1 Rabbit Polyclonal Antibody

CAB15216



Product Information

Product SKU:	CAB15216	Gene ID:	202500	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	66kDa	Conjugate:	Unconjugated
Calculated MW:	56kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in flagellated sperm motility. Predicted to be located in sperm flagellum.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	DRC5; D6S46; FAP155; TCTE1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TCTE1 (NP_872345.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.