

Product Information

Product SKU:	CAB15225	Gene ID:	344805	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	94kDa	Conjugate:	Unconjugated
Calculated MW:	94kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable serine-type peptidase activity. Predicted to be involved in proteolysis. Predicted to be integral component of membrane. Predicted to be active in plasma membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	TMPRSS7
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 478-717 of human TMPRSS7 (NP_001036040.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.