

ligp1 Rabbit Polyclonal Antibody



CAB15226

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

47kDa

Calculated MW:

47kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

60440

Uniprot

Synonyms:

2900074L10Rik; AI046432; AW111922; Ifgga1; ligp; Irga6; ligp1


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 150-250 of mouse ligp1 (NP_068564.4).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - ligp1 Rabbit pAb (CAB15226)