

Product Information

Product SKU:	CAB15228	Gene ID:	79674	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	95kDa	Conjugate:	Unconjugated
Calculated MW:	95kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable phosphatidylinositol-5-phosphate binding activity. Involved in negative regulation of SMAD protein signal transduction and negative regulation of transforming growth factor beta receptor signaling pathway. Predicted to be active in plasma membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	MELT; VEPH; VEPH1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 448-547 of human VEPH1 (NP_078897.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.