
Product Information

Product SKU:	CAB15241	Gene ID:	285671	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	65kDa	Conjugate:	Unconjugated
Calculated MW:	68kDa	Isotype:	IgG

Immunogen Information

Background:	Predicted to enable ubiquitin conjugating enzyme binding activity and ubiquitin protein ligase activity. Predicted to be involved in norepinephrine metabolic process; positive regulation of proteasomal ubiquitin-dependent protein catabolic process; and serotonin metabolic process. Predicted to act upstream of or within several processes, including adult behavior; positive regulation of protein ubiquitination; and protein polyubiquitination. Predicted to be located in nuclear envelope. Predicted to be integral component of membrane. Predicted to be intrinsic component of endoplasmic reticulum membrane.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	RINES; RNF180
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-416 of human RNF180 (NP_848627.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.