

SIK1 Rabbit Polyclonal Antibody



CAB15247

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

84kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

150094

Uniprot

P57059

Synonyms:

SIK1; MSK; SIK; SIK-1; SNF1LK; EIEE30

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SIK1 (NP_775490.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
