

Product Information

Product SKU:	CAB15291	Gene ID:	4618	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	27kDa/37kDa	Conjugate:	Unconjugated
Calculated MW:	27kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a probable basic helix-loop-helix (bHLH) DNA binding protein involved in muscle differentiation. The encoded protein likely acts as a heterodimer with another bHLH protein. Defects in this gene are a cause of autosomal dominant centronuclear myopathy (ADCNM).

Recommended Dilution:	WB,1:200 - 1:2000
Synonyms:	CNM3; MRF4; myf-6; bHLHc4; MYF6
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 80-180 of human MYF6 (NP_002460.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.