

Product Information

Product SKU:	CAB15293	Gene ID:	4622	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	250kDa	Conjugate:	Unconjugated
Calculated MW:	223kDa	Isotype:	IgG

Immunogen Information

Background: Enables double-stranded RNA binding activity. Involved in muscle contraction. Located in myofibril.

Recommended Dilution:	WB,1:100 - 1:500 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200
Synonyms:	MYH2B; MyHC-2B; MyHC-IIb; MYH4
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1500-1600 of human MYH4 (NP_060003.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.