

CAB15319

Product Information

Product SKU:	CAB15319	Gene ID:	6769	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	44kDa	Conjugate:	Unconjugated
Calculated MW:	45kDa	Isotype:	IgG

Immunogen Information

Background:	Predicted to enable transmembrane transporter binding activity. Predicted to be involved in positive regulation of protein localization to plasma membrane; positive regulation of voltage-gated calcium channel activity; and skeletal muscle contraction. Predicted to act upstream of or within cellular response to heat; muscle contraction; and regulation of voltage-gated calcium channel activity. Predicted to be located in T-tubule. Predicted to be extrinsic component of cytoplasmic side of plasma membrane.
Recommended Dilution:	WB,1:200 - 1:2000
Synonyms:	STAC1; STAC
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 153-402 of human STAC (NP_003140.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.