

TPTE Rabbit Polyclonal Antibody



CAB15325

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

49kDa/60kDa/62kDa/64kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a PTEN-related tyrosine phosphatase which may play a role in the signal transduction pathways of the endocrine or spermatogenic function of the testis. Alternative splicing results in multiple transcript variants.

Immunogen information

Gene ID:

7179

Uniprot

P56180

Synonyms:

TPTE; CT44; PTEN2

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 459-533 of human TPTE (NP_954868.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
