

DGAT1 Rabbit Polyclonal Antibody



CAB15340

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

55kDa

Applications:

WB

Reactivity:

Mouse, Rat

Protein Background

This gene encodes an multipass transmembrane protein that functions as a key metabolic enzyme. The encoded protein catalyzes the conversion of diacylglycerol and fatty acyl CoA to triacylglycerol. This enzyme can also transfer acyl CoA to retinol. Activity of this protein may be associated with obesity and other metabolic diseases.

Immunogen information

Gene ID:

8694

Uniprot

O75907

Synonyms:

DGAT1; ARAT; ARGP1; DGAT; DIAR7

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human DGAT1 (NP_036211.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
