

FNDC3A Rabbit Polyclonal Antibody



CAB15404

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

125kDa/131kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

22862

Uniprot

Q9Y2H6

Synonyms:

FNDC3A; FNDC3; HUGO; bA203I16.1; bA203I16.5

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 900-1170 of human FNDC3A (NP_001265367.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
