

MPRIP Rabbit Polyclonal Antibody



CAB15409

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

114kDa/116kDa/118kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23164

Uniprot

Q6WCQ1

Synonyms:

MPRIP; M-RIP; MRIP; RHOIP3; RIP3; p116Rip

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human MPRIP (NP_055949.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
