

CAB15419

Product Information

Product SKU:	CAB15419	Gene ID:	25926	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	80-85kDa	Conjugate:	Unconjugated
Calculated MW:	81kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables RNA binding activity. Involved in maturation of SSU-rRNA and positive regulation of transcription of nucleolar large rRNA by RNA polymerase I. Located in nucleolus. Part of t-UTP complex.

Recommended Dilution:	WB,1:200 - 1:2000
Synonyms:	NOL11
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 300-600 of human NOL11 (NP_056277.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.