

ZNF451 Rabbit Polyclonal Antibody



CAB15422

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

63kDa/115kDa/121kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

26036

Uniprot

Q9Y4E5

Synonyms:

ZNF451; COASTER; dJ41711.1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 600-900 of human ZNF451 (NP_001026794.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
