

---

## Product Information

<b>Product SKU:</b>	CAB15422	<b>Gene ID:</b>	26036	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	Refer to figures	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	121kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Enables SUMO ligase activity and transcription corepressor activity. Involved in negative regulation of nitrogen compound metabolic process; negative regulation of transforming growth factor beta receptor signaling pathway; and protein sumoylation. Located in nucleus.

**Recommended Dilution:** WB,1:200 - 1:2000

**Synonyms:** ZATT; COASTER; dJ41711.1; ZNF451

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 600-900 of human ZNF451 (NP\_001026794.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.