

THEG Rabbit Polyclonal Antibody



CAB15447

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

40kDa/43kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

This gene is specifically expressed in the nucleus of haploid male germ cells. The orthologous gene in mice encodes a protein that may play a role in protein assembly through interactions with T-complex protein 1 subunit epsilon. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Immunogen information

Gene ID:

51298

Uniprot

Q9P2T0

Synonyms:

THEG; CT56; THEG1

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-379 of human THEG (NP_057669.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
