

SPIN2A Rabbit Polyclonal Antibody



CAB15456

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

29kDa

Applications:

WB

Reactivity:

Human, Mouse

Protein Background

This gene encodes one of three members of the DXF34 gene family, located in a 100-kb region of chromosome Xp11.21.

Immunogen information

Gene ID:

54466

Uniprot

Q99865

Synonyms:

SPIN2A; DXF34; SPIN2; TDRD25; dJ323P24.1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-50 of human SPIN2A (NP_061876.3).

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
