

Product Information

Product SKU:	CAB15482	Gene ID:	57092	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	26kDa	Conjugate:	Unconjugated
Calculated MW:	19kDa	Isotype:	IgG

Immunogen Information

Background: Involved in proteasome-mediated ubiquitin-dependent protein catabolic process and protein ubiquitination. Located in nuclear body.

Recommended Dilution:	WB,1:200 - 1:2000
Synonyms:	PCNP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human PCNP (NP_065090.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.