

KIAA1324 Rabbit Polyclonal Antibody



CAB15490

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

101kDa/109kDa/110kDa/111kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Expression of this gene is induced by estrogen and the encoded protein has been characterized as a transmembrane protein. The encoded protein has been found in to correlate with survival in certain carcinomas (PMID: 21102415) and may be important for cellular response to stress (PMID: 21072319). Alternative splicing results in multiple transcript variants.

Immunogen information

Gene ID:

57535

Uniprot

Q6UXG2

Synonyms:

KIAA1324; EIG121

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 42-300 of human KIAA1324 (NP_001253977.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
